

# technical memorandum

## Daresbury Laboratory

DL/SCI/TM38E

X-RAY ABSORPTION SPECTROSCOPY

by

G.P. DIAKUN, G.N. GREAVES, S.S. HASNAIN and P.D. QUINN, Daresbury Laboratory

JANUARY, 1984

Science &amp; Engineering Research Council

Daresbury Laboratory

Daresbury, Warrington WA4 4AD

LENDING COPY

© SCIENCE AND ENGINEERING RESEARCH COUNCIL 1984

Enquiries about copyright and reproduction should be addressed to:—  
The Librarian, Daresbury Laboratory, Daresbury, Warrington,  
WA4 4AD.

**IMPORTANT**

The SERC does not accept any responsibility for loss or damage arising from the use of information contained in any of its reports or in any communication about its tests or investigations.

	<u>Page</u>
INTRODUCTION.	2
Chapter 1 - STATION SOFTWARE.	5
1.1 Initial check.	5
1.2 Loading the various programmes.	5
1.3 Transferring data to the mainframe.	8
1.4 Running procedure for transmission measurements.	9
1.5 Running procedure for fluorescence measurements.	19
Chapter 2 - STATION HARDWARE	23
2.1 Setting up and operation of order sorting monochromator.	23
2.2 Setting up and operation of fluorescence detectors.	27
2.3 Operation with toroidal mirror.	29
2.4 Guidelines for EXAFS measurements.	30
Chapter 3 - DATA REDUCTION AND ANALYSIS	35
3.1 EXCALIB	35
3.2 EXBACK	37
3.3 EXCURVE	40
APPENDIX	43

This manual describes facilities for x-ray absorption spectroscopy at the SRS in the spectral range 0.3 - 3.2 Å or 40 - 4 keV. It is divided into three main chapters. The first covers the station software (1.1 to 1.5), the second the station hardware (2.1 to 2.4) and the third is devoted to the analysis programmes.

In particular Chapter 1 describes loading the setting-up and data acquisition programmes (1.2) and transferring data from the experimental station to the mainframe (1.3). A command-by-command description of the software is then given - both for transmission as well as for fluorescence measurements (1.4 and 1.5). Chapter 2 covers the operation of the major items of equipment on the station: the order-sorting monochromator (2.1) the scintillation PMTs for fluorescence detection (2.2) and the Pt-coated toroidal mirror (2.3). Chapter 2 also includes broad guidelines for making both transmission and fluorescence EXAFS measurements (2.4). Finally, in Chapter 3 an introduction is provided to the data reduction and EXAFS analysis programmes available in the SRS Programme Library. The three principal programmes EXCALIB (3.1), EXBACK (3.2) and EXCURVE (3.3) are described.

A schematic of the EXAFS equipment is shown in figure 1. Unfocused or focused white x-rays from the beryllium window pass through monochromator vessels where they are defined in a vertical angle and horizontal aperture and monochromated. The EXAFS monochromator operates under rough vacuum pressure. Monochromated x-rays leave the exit slit vessel through a further beryllium window. EXAFS measurements are usually made in air either by monitoring

- (a) transmission through the sample
- or (b) fluorescence emitted from the sample.

The user can expect that the entrance and exit slits will be aligned with respect to the monochromator and that the Bragg angle of the monochromator will be approximately calibrated. Ion chambers will be in position and have suitable rare gas mixtures for the particular wavelength the monochromator has been left at. If fluorescence detectors are required help should be sought from the station master. Precise calibration of the mono-

chromator is left to the user and requires the use of a metal foil whose K and L edges fall in the wavelength range of interest. Monochromatic beam can be detected using Polaroid film and a laser is provided to enable the sample to be aligned to the monochromated beam. Various cryostats and a furnace are available. If these or other environmental stages are required the Station Master or Local Contact should be advised in advance of the scheduled beam time.

New users should be aware that x-ray beam lines at the SRS are multi-station ports and that there are two kinds of x-ray shutter. The port shutter, which controls the light to all the stations is under the control of the crew in the Main Control Room (MCR). Should this close for any reason the appropriate light will go out and the MCR should be contacted for it to be re-opened (telephone 560, 561 or 562). Persistent port shutter problems must be reported to the Beam Port Coordinator as soon as possible. The station shutter is under the control of the experimenter. Once the hutch is correctly interlocked the station shutter can be opened and closed as desired. See the instruction notice on the hutch for details of the hutch interlocking procedure. Station shutter problems should be reported to the Station Master.

The documentation that follows is designed to be reasonably comprehensive and will be updated as facilities at the SRS develop. If problems arise during data collection or analysis that are outside the scope of this manual, the Station Master or the Local Contact should be advised. In any event it is extremely important that all difficulties and problems encountered by users should be fully recorded in the station log-book at the time.

PEOPLE TO CONTACT FOR HELP

		Ext. No.	Bleep
S.S. Hasnain	(Station Master 7.1)	273	225
G.P. Diakun	(Deputy Station Master 7.1)	273	
G.P. Diakun	(Station Master 9.4)	273	
G.N. Greaves	(Deputy Station Master 9.4)	335	107
P.A. Ridley	CSE	405	105
E. Pantos	CSE, SRS Programme Library	275	
G. Mant	CSE	275	

J. Sheldon	Mechanical Services (e.g. gas bottles, plumbing, rotary pumps)	557	173
		Ext. No.	Bleep
User Support Team		560,561,562	
Beam Port Coordinators:			
X-ray 7	G.P. Diakun (602124) (Deputy G.N. Greaves)		
X-ray 8	S.S. Haanain (93-710843) (Deputy J.S. Morgan)		
Wiggler 9	G.N. Greaves (935-63194) (Deputy P.J. Duke)		

1.1 INITIAL CHECKS

1.1.1) Check the CAMAC power supply is switched on.

1.1.2) Check the disk drivers are switched on.

1.1.3) Programme status. Press 'carriage return'. If RT-11 is loaded the VDU will respond with the prompt '.', if CATEX is also loaded the prompt '+' will appear on the screen, if the EXAFS programme is loaded

INVALID COMMAND FOR HELP TYPE H  
COMMAND?

will appear on the console.

1.2 LOADING THE VARIOUS PROGRAMMES

The EXAFS and system programmes are stored at the station on a Winchester Disk which also collects a permanent record of the experimental data collected. The loading procedure depends on the state the system is in. RT-11 responds with a '.' prompt, CATEX with a '+' prompt and the EXAFS programme with a 'COMMAND?' prompt. The loading procedure should therefore be taken up from the prompts shown on the screen when carriage return is pressed. The full procedure is necessary if the system crashes for some reason (e.g. power failure).

1.2.1) Using the keys on the VDU

Press BREAK  
Enter 762000G

The system run light should now be illuminated, if the system does not load or the run light goes off repeat the sequence. If this fails seek advice.

1.2.2) When the system has been loaded the console will respond with:

RT11 FB(S) V04.00  
DATE  
TIME

When the RT-11 prompt '.' has appeared CATEX may be loaded by typing  
CAT 'carriage return'

The system responds with:

\*\*\*\*\*CATEX Version 2\*\*\*\*\*  
LAST MODIFICATION <date>  
+

When the CATEX prompt '+' has appeared the data acquisition programme may be loaded by typing either

IRD SY:LOAD 'carriage return'

for transmission measurements or

IRD SY:LOAD2 'carriage return'

for fluorescence measurements.

If at any time during the running of the programme an unforeseen error occurs and the CATEX prompt appears the system should be reloaded with the above commands. Important: please report all errors in EXAFS7 logbook. When either programme has been loaded the console responds with:

TABLE HEIGHT (MICRONS)=xxxx (e.g. 16000)  
MONOCHROMATOR POSITION=xxxx (e.g. 24000)

The table height is the height of the optical bench (containing exit slits, ion chambers etc) that is offset from straight through beam. This offset is required because the white radiation when monochromatised is displaced by  $2D\cos\theta$ , where D is the separation between the two crystal surfaces of the monochromator and  $\theta$  is the angle of the monochromator with respect to the incident beam, or the Bragg angle. The monochromator position is given in thousandths of a degree, the value 24000, therefore, is equivalent to 24 degrees. When moving the monochromator one step this is equal to a thousandth of a degree.

A list of commands also appears on the console. For the programme LOAD these are:

TYPE L LIST SCAN PARAMETERS

C CHANCE SCAN PARAMETERS  
E MOVE SLITS, TABLE, EXPERIMENT AND MIRROR  
P MONOCHROMATOR POSITION  
M MOVE MONOCHROMATOR  
S START DATA COLLECTION  
I FOR INFORMATION ON K AND L(3) EDGES  
R READ ION CHAMBERS  
T HAVE TABLE SYNCHRONISED WITH MONOCHROMATOR  
V TABLE HEIGHT  
O MOVE COARSE MOTOR  
D % HARMONIC REJECTION FOR DC  
H HELP  
COMMAND?

For the programme LOAD2 these are:

TYPE L LIST SCAN PARAMETERS

C CHANCE SCAN PARAMETERS  
E MOVE SLITS, TABLE, EXPERIMENT AND MIRROR  
P MONOCHROMATOR POSITION  
M MOVE MONOCHROMATOR  
S START DATA COLLECTION  
I FOR INFORMATION ON K AND L(3) EDGES  
R READ ION CHAMBERS  
T HAVE TABLE SYNCHRONISED WITH MONOCHROMATOR  
V TABLE HEIGHT  
O MOVE COARSE MOTOR  
D % HARMONIC REJECTION FOR DC  
H HELP  
A TO SET UP FLUORESCENCE DETECTORS  
COMMAND?

These commands enable a user to set up his experimental parameters prior to data collection. They allow adjustments to be made to the slits, optical bench, mirror and the Bragg angle of the monochromator. They enable the counts in the reference and signal ion chambers to be read and the energy window of the fluorescence detectors to be set. They also provide the facility for the user to set up spectral regions for a scan. The commands are described in detail in the next section.

1.2.3) Important: ensures that the Station Master has updated the system with your current grant number.

1.3 TRANSFERRING DATA TO THE MAINFRAME

Once a spectrum is complete, data can be transferred from the experimental station to the SRS mainframe. This is useful as it enables initial data to be analysed as the experiments progress. In order to transfer the data type 'F' to finish session in the scanning programme. The command appears in the transfer to the scanning programme by first using command 'S'. Then type

!AB ' carriage return'

The console will respond with the prompt '.' which means the user is in RT77.

Once in RT11, load DQ and SL, viz

.LOAD DQ

.LOAD SL

CATEX should then be loaded by typing

CAT 'carriage return'

The VDU will respond with

\*\*\*\*\*CATEX Version 2\*\*\*\*\*

LAST MODIFICATION <data>

+

As soon as the CATEX prompt '+' has appeared type

!RS SY:CATSTR

A prompt '+' will appear when the subroutine has been loaded, then enter

!RD SY:TFILE

This loads the programme which allows data to be transferred to the NAS7000. The CATEX prompt '+' is displayed on the VDU after the programme has been read in. In order to obtain a list of the files on the data disk so that the user can choose which to transfer, type

!DI

and the directory of the data disk is then printed on the VDU.

To start the programme type

15T

and then reply to the questions as they appear. When files have been transferred unload DQ and SL, viz

+IAB

.UNLOAD DQ

.UNLOAD SL

Then load CATEX and read the LOAD programme as described above in section 1.2.

#### How to Re-initialise a Floppy Disk

If for any reason Floppy Disks are used for data collection these are re-usable. In order to re-initialise a Floppy Disk transfer to RT-11. To enter RT-11 from CATEX (which has the prompt '+') type

IAB

The VDU will respond with the prompt '.'. To wipe the disk clean type

INI/BA FD1: 'carriage return'

The system will ask if you are sure and the reply is

Y 'carriage return'

The system will then clean the disk and check if any bad blocks are present. If this is so do not re-use the disk to collect data. Also do not type 'INI/BA FDO:' to re-initialise as this will wipe the system disk clean. To get back to CATEX after re-initialising type

CAT 'carriage return'

#### 1.4 RUNNING PROCEDURE FOR TRANSMISSION MEASUREMENTS

A transmission measurement involves placing the sample between the two ion chambers and scanning the monochromator through the x-ray absorption edge of one of the elements present. Figure 2 shows the transmission spectrum of Ni foil. The specimen should be of appropriate thickness (2.4.2) and the ion chambers should be filled with rare gas mixtures to suit the spectral wavelength range chosen (2.4.8). Details of K and L edges - wave-

lengths, energiss, monochromator Bragg angles and suggested spectral ranges are provided in the INFORMATION command 'I' described below.

With the programme LOAD installed, typing 'H' displays the commands necessary for setting up a transmission EXAFS spectrum.

#### TYPE L LIST SCAN PARAMETERS

C CHANCE SCAN PARAMETERS

E MOVE SLITS, TABLE, EXPERIMENT AND MIRROR

P MONOCHROMATOR POSITION

M MOVE MONOCHROMATOR

S START DATA COLLECTION

I FOR INFORMATION ON K AND L(3) EDGES

R READ ION CHAMBERS

T HAVE TABLE SYNCHRONISED WITH MONOCHROMATOR

V TABLE HEIGHT

O MOVE COARSE MOTOR

D HARMONIC REJECTION FOR DC

H HELP

COMMAND?

The user may enter any of the commands available by typing the appropriate letter (no 'carriage return necessary). If a letter is typed which is not on the list of commands the programme responds with:

INVALID COMMAND FOR HELP TYPE H

COMMAND?

#### COMMANDS

##### 1.4.1) L - LIST SCAN PARAMETERS

This lists the number of scans and regions (for each scan) previously entered into the programme and these may be changed using command 'C'. The contents of each region are also displayed on the VDU, for example:

SCANS 1	REGIONS 2
REGION 1	
START 14100	TIME 100000
INCR 10	REF 8000000
END 13000	SIG 80000
REGION 2	
START 13000	TIME 100000

```
INCR      2      REF 8000000
END 12000      SIG 100000
COMMAND?
```

These scan parameters are for a single scan of two regions starting at 14.1 degrees and ending at 12.0 degrees. The number of steps the monochromator moves before a reading is taken is determined by inor (short for increment). Time is measured in milliseconds and the values for the reference and signal are counts from the ion chambers. The usual mode of setting up the regions is to give large values for time and reference so that the counting is limited on signal, i.e. once the number for the signal is reached the accumulation of counts for references and time are stopped simultaneously and a reading taken, the monochromator then moves to the next angle. In this way the reference counts are inversely proportional to the sample transmission. Command 'R' can be used to choose appropriate values for time, reference and signal.

#### 1.4.2) C - CHANGE SCAN PARAMETERS

Typing this letter allows a user to change the number of scans and regions he wishes to have from those already inputted into the programme. The parameters for each region can also be changed. The VDU responds with:

```
TYPE P CHANGE PARAMETERS
R CHANGE REGION
E EXIT
```

By typing 'P' the console prompts the user for the number of scans he wishes to run repetitively, the maximum number allowed is 100. If a value greater than this is entered the programme will use the default of 100. If the user does not wish to alter the number already present then 'carriage return' is pressed. The user is then asked to input the number of regions he requires for a scan, the maximum number allowed is 10. If a value greater than this is inputted the programme assumes the default of 10. After the value has been entered the user is taken back to the above three commands.

Typing 'R' causes the user to be prompted for a region number. If he enters a value larger than the regions specified, the programme responds with:

```
INVALID REGION
?
```

By entering an appropriate region number the console displays:

```
B START OF REGION
I INCREMENT
F END OF REGION
T CHANGE TIME
R REF.I.C.
S SIG.I.C.
N NEW REGION
E EXIT
```

R, I and F are used for entering the start, increment and end of a region in a spectrum in millidegrees. The start of the region should have a larger number than the end of a region. T, R and S are used as limits on time (milliseconds), reference ion chamber (in counts) and signal ion chamber (in counts) respectively. The maximum allowed time is 100,000, and for reference and signal ion chamber the maximum counts are 8,300,000. When a user exits from these sub-commands the programme returns to the main commands. An example is given below of a three region scan for the copper edge of a copper foil using a Si(111) monochromator.

```
REGION 1
START 13400      TIME 100000
INCR      20      REF 8000000
END 12740      SIG 500000
REGION 2
START 12740      TIME 100000
INCR      1      REF 8000000
END 12490      SIG 100000
REGION 3
START 12490      TIME 100000
INCR      3      REF 8000000
END 11800      SIG 100000
```

The first region has a coarse increment and finishes about 20 steps before the start of the copper edge, the second region encompasses the edge and near edge structure, and has the finest increment. The final region



covers the EXAFS oscillations and extends to about 800 eV above the edge. Because the sample absorbs in regions 2 and 3 the number of counts entered for the signal is lower than for region 1. This takes into account the decrease in x-rays entering the signal ion chamber above the absorption edge so that the counting time for each point remains roughly the same for all regions. As the EXAFS oscillations are broader than the near edge structure the increment in region 2 is smaller than in region 3.

#### 1.4.3) E - MOVE SLITS, TABLE, EXPERIMENT AND MIRRORS

Typing this letter the console responds with:

A MOVE ENTRANCE SLIT VERTICALLY  
B MOVE ENTRANCE SLIT HORIZONTALLY  
C MOVE EXIT SLIT VERTICALLY  
D MOVE EXIT SLIT HORIZONTALLY  
F MOVE FRONT OF TABLE  
G MOVE BACK OF TABLE  
H MOVE FRONT OF EXPERIMENT  
I MOVE BACK OF EXPERIMENT  
J MOVE FRONT OF MIRROR  
K MOVE BACK OF MIRROR  
E EXIT

This command allows the user to perform the following adjustments:

- 1) Open and close the slits before and after the monochromator by entering a positive or negative number (A-D).
- 2) Raise and lower the optical bench to compensate for any change in position of monochromatic beam (F-G). When the table is moved its new height is displayed in microns. The number is calibrated with respect to the white radiation passing through the entrance and exit slits with no deviation encountered. A positive value for the table height means the optical bench is displaced upwards as compared to the straight through beam. The front and back of the table can be moved independently to facilitate levelling of the equipment. The value for the table height is taken from the front of the table, so remember to move the back of the table to retain the bench in a horizontal position.

- 3) Raise and lower the whole of the experimental assembly in the EXAFS hutch to compensate for the reflected beam off the mirror when placed in the synchrotron radiation (H-I). This adjustment should not be made without permission of the station master.
- 4) The mirror to be moved in and out of the x-ray beam and also variation in the glancing angle of the mirror (section 2.3). This adjustment should not be made without permission of the station master.

#### 1.4.4) P - MONOCHROMATOR POSITION

This informs the user of the present position of the monochromator and allows him to recalibrate if required, if not 'carriage return' is pressed. Calibration can not be made by driving the monochromator to the 'turning point' of the absorption edge of the element and entering the angle tabulated in command 'I'.

#### 1.4.5) M - MOVE MONOCHROMATOR

This command enables the monochromator to be increased or decreased in Bragg angle by a desired number of steps (1 step is 1 millidegree). The maximum number of steps the monochromator may have at any one time is 30000, if a value greater than 30000 is entered the console informs the user of the maximum permissible steps he is allowed and prompts for a new value to be typed. When the monochromator has finished moving the VDU responds with its new angle and also the table height.

If the wrong direction is entered by mistake do not interrupt until the movement has been completed and re-enter in the opposite direction with double the number of steps to correct. Zero steps is not allowed, the console will prompt the user for a new value.

#### 1.4.6) I - FOR INFORMATION ON K AND L(3) EDGES

When 'I' is inputted the programme prompts the user for the particular edge he is interested in:

ENTER 1 FOR K EDGE 2 FOR L(3) EDGE

After the appropriate value is entered the console responds with:  
ENTER ELEMENT

All that is required is for the chemical symbol of the element concerned to be entered. Should the chemical symbol have only one letter the user should type this letter followed by a blank then 'carriage return', if it has two letters no space is required just 'carriage return'. The programme responds with information on the edge in terms of where it occurs (in energy (KeV), wavelength (Å) and monochromator position for Si(220) and Si(111) crystals. Sample scan ranges of 300 eV before the edge and 600 eV after the edge expressed as monochromator steps are also given.

#### 1.4.7) R - READ ION CHAMBERS

With this command a user may read the counts in the ion chambers for a fixed time during the setting up procedure. The default value for time is 1000 milliseconds and this may be altered by the user. The reference and signal counts both have default values of 1,000,000 which are fixed. This system is useful for setting up the reference and signal limits in the scan parameters. With the shutters closed the user is also able to obtain dark current (amplifier off-set) readings.

#### 1.4.8) T - HAVE TABLE SYNCHRONISED WITH MONOCHROMATOR

Because the separation between the two crystals is fixed (D) the output beam from the monochromator moves relative to the entrance beam as  $\theta$  is scanned, the displacement is given by  $2D\cos\theta$ . The 'T' option allows the user to track the output beam with the exit slits and optical bench. The user is prompted to input the crystal separation (millimeters) and the amount of beam movement permissible before the table is moved to correct for it (microns). The programme responds by converting the beam movement into steps for the stepper motors on the optical bench.

WHAT SEPARATION BETWEEN CRYSTAL SURFACES (MM) <VALUE>

Input new value 'carriage return'

AMOUNT OF BEAM MOVEMENT (MICRONS) BEFORE TABLE MOVED <VALUE>

Input new value 'carriage return'

STEP=?

If a user does not wish to track the beam then zero should be entered for both separation of crystal surfaces and amount of beam movement. The separation of the crystal surfaces of the two crystal order-sorting monochromator is 15 mm and the separation for the Si(111) channel cut (D-shaped) is 4 mm for other crystals, ask the station master. A reasonable step size

for the synchronised motion is 1-5 microns.

#### 1.4.9) V - TABLE HEIGHT

This command informs the user of the present height of the optical bench with respect to the straight through beam and may be recalibrated if required. If the monochromator angle is calibrated, the table height is given by  $2D\cos\theta$ , where D is the separation of the crystal surfaces.

#### 1.4.10) O - MOVE COARSE MOTOR

This command should only be used when a double crystal monochromator is being operated. It allows a stepper motor to turn which alters the angle of the first crystal with respect to the second crystal. This varies the intensity of monochromatic light leaving the monochromator and with it the degree of harmonic contamination. Fine control is effected with the solenoid drive which can be operated manually (see '2.1 SETTING UP AND OPERATION OF ORDER SORTING MONOCHROMATOR').

#### 1.4.11) D - % HARMONIC REJECTION FOR D.C. (Double Crystal)

This command should only be used if the double crystal monochromator is installed. It allows a user to carry out a scan with approximately the same harmonic rejection throughout. This is done by a servo linkage between the output of the reference ion chamber (I(O)) and the angle of the first crystal of the monochromator. Before using this command the peak height for I(O) at the beginning and end of the scan should be noted. These can be obtained using the coarse motor (command 'O') in conjunction with the solenoid manual drive control (see '2.1 SETTING UP AND OPERATION OF ORDER SORTING MONOCHROMATOR').

When 'D' is typed the console inquires the percentage of I(O) desired for the scan. This is directly related to the degree of harmonic rejection. For Si(220) on EXAFS station 7.1

Scanning around the 3 Å region the % of I(O) should be 40%

Scanning around the 2 Å region the % of I(O) should be 50%

Scanning around the 1 Å region the % of I(O) should be 75%

The above values are meant as rough guidelines for the SRS ring operating at 2 GeV. If the storage ring is running at a lower energy these may

be increased. However values of greater than 85% I(O) should not be used as the monochromator servo mechanism will have difficulty holding a level so close to the rocking curve peak during the scan. As an example fig.3 shows the effects of harmonic contamination on an absorption spectrum (in this case V foil).

After entering the % of I(O) required the programme prompts the user for the rocking curve peak height of I(O) at the start of the scan, followed by the rocking curve peak height at the end of the scan range. When these values have been entered the programme returns to the list of main commands. Also the user will observe a value displayed on the servo reference on the right-hand DVM of the EXAFS servo panel. This displays the I(O) value which will act as the servo reference at the start of the scan. In order to limit the solenoid current, the monochromator should be detuned by using the coarse motor motion (command 'O'). The output of the reference ion chamber displayed on the left-hand DVM on the 'EXAFS CURRENT MONITOR' panel should be adjusted as close as possible to the servo reference value shown above. The servo current is displayed on the left-hand DVM of the EXAFS servo panel (see '2.1 SETTING UP AND OPERATION OF ORDER SORTING MONOCHROMATOR').

#### 1.4.12) S - START DATA COLLECTION

This allows a user to start collecting an EXAFS spectrum with the data stored on the winchester disk. However before issuing this command:

- 1) Check that the position of the monochromator is greater or equal to the start of the first region otherwise the scan will automatically abort.
- 2) Ensure the end of one region does not overlap with the start of the following region (they may, however, have the same value).
- 3) Check the start of a region is always greater than the end of the same region.

Issuing the command 'S' causes the parameters to be saved on the winchester disk and the setting up programme is then replaced by the scanning programme. This procedure takes approximately 5 seconds to perform, so wait till the console responds with the following commands before pressing any keys on the VDU.

```
TYPE S START DATA COLLECTION
C CHANGE SETTING UP PROCEDURE
F FINISH SESSION
```

When the letter 'S' is typed the console responds with  
EVERYTHING OK Y/N

If the user is satisfied he is then prompted for an expansion factor for the y-axis of the scan which is to be plotted on the session VDU.

```
EXPANSION FOR Y AXIS 3.0 LARGER INTEGER SMALLER EXPANSION
?
```

For a metal foil the value entered is typically between 15 to 40 and for a dilute sample 1 to 3. This is a cosmetic to the display and does not affect the data collected. After the scaling factor has been entered the console displays

```
STARTING RUN NUMBER R*****
NUMBER OF BLOCKS ALLOCATED TO THE DATA FILE = ****
ENTER TITLE
ENTER CONDITION
ENTER COMMENT
```

After the title has been entered (maximum of 40 characters) the user is prompted for conditions 1 to 3 (up to 8 characters may be entered for each condition), and is then prompted for comments (up to 40 characters). To exit from this sequence press 'carriage return'.

The dataset produced on the winchester disk will have its run number with the letter 'R' before it and will have the title R\*\*\*\*.DAT on the disk. The dataset will contain the standard SRS header before the data, e.g.:-

```
€SRS
SRSRUN=1677,SRSDAT=820421,SRSTIM=171441
SRSSSTN='EXP1',SRSPRJ='SR1234',SRSEXP='GLASSES'
SRSTLE='*****'
CONDITION 1='*****'
COMMENT='*****'
€END
```

SRSRUN, SRSDAT and SRSTIM informs the user of the run number of his data, the year, month and day when the data was collected and the time at which data collection started. SRSSSTN is the station where the experiment was performed. SRSPRJ and SRSEXP are respectively the project number and

the type of system under investigation. After this header the title with any conditions and comments are printed followed by the data.

If you wish to abort during a scan, type 'A' followed by 'carriage return'.

The data collected so far will be saved on winchester disk and at the end of the file the user will be informed that the scan was aborted.

During a multiple scan, after entering the relevant details for the initial scan, the user is not required to input the same information for subsequent scans. At the end of a scan, or multiple scans, the user has the option to compare the last three datasets collected to be plotted on the session screen.

Typing 'C' causes the setting up programme to replace the scanning programme, the user can then make alterations to his input parameters using command 'C'.

Typing 'P' removes the compiled programme, however the last input parameters the user used are saved on winchester disk. To re-load the programme after letter 'P' has been typed, input

```
!RD SY:LOAD
```

(see '1.2 LOADING THE VARIOUS PROGRAMMES').

### 1.5 RUNNING PROCEDURE FOR FLUORESCENCE MEASUREMENTS

In order to make a fluorescence measurement the sample should be placed behind the reference ion chamber, inclined at 45 degrees with respect to the monochromatic beam and surrounded by an array of scintillation detectors. A photograph of the detector arrangement used for fluorescence measurements is shown in fig.4. An x-ray excitation spectrum is obtained by scanning the monochromator through the absorption edge of the dilute component. The specimen should be of appropriate thickness (2.4.2). The fluorescence detector should be arranged (with filters if necessary) to optimise the fluorescence to scatter (see section 2.2). Fluorescence spectra

for a selection of concentrations of  $\text{CuSO}_4$ , are presented in fig.5.

Having loaded the LOAD2 programme typing 'H' displays the commands necessary for setting up a fluorescence EXAFS spectrum

```
TYPE L LIST SCAN PARAMETERS
C CHANGE SCAN PARAMETERS
E MOVE SLITS, TABLE, EXPERIMENT AND MIRROR
P MONOCHROMATOR POSITION
M MOVE MONOCHROMATOR
S START DATA COLLECTION
I FOR INFORMATION ON K AND L(3) EDGES
R READ ION CHAMBERS
T HAVE TABLE SYNCHRONISED WITH MONOCHROMATOR
V TABLE HEIGHT
O MOVE COARSE MOTOR
D * HARMONIC REJECTION FOR DC
H HELP
A TO SET UP THE FLUORESCENCE DETECTORS
COMMAND?
```

The majority of these commands are identical to those described previously for transmission measurements (1.4). The user may enter any of the commands by typing the appropriate letter - without a carriage return. If a letter is typed which is not included in the list above the programme responds with:

```
INVALID COMMAND FOR HELP TYPE H
COMMAND?
```

#### COMMANDS

##### 1.5.1) A - TO SET UP FLUORESCENCE DETECTORS

The idea of this small routine is to allow a user to set up the detectors individually and then to match them with each other. Also it allows a user to choose his MCA window for data collection. Note that at the start of each day the detectors should be checked individually to ensure they have not drifted or broken down. Below more detailed information is given about the routine.

Typing 'A' removes the main programme from memory and also stores any

changes on the winchester disk. The programme asks the user

'DO YOU WISH TO RETURN TO MAIN PROGRAMME'

If 'N' is typed the VDU responds with

'TIME'

The console waits for a value to be inputed. (Important a value of 50 = 1 second).

After counting for the time specified the programme asks

'DO YOU WISH TO COMPARE SPECTRA' Y/N

If the answer is 'Y' the terminal displays

'J=x'

x can be 1 or 2. This corresponds to a memory allocation for displaying two spectra on the Hytec TV screen which allows you to set up the detectors to match one another. For instance using memory 1 for the correctly adjusted detector, the other detectors can be matched to this via memory 2.

If the reply is 'N' the spectrum is placed in memory 1 (J=1) and is displayed on the Hytec TV screen with a line of dots at regular intervals half way up the screen. These dots are spaced at every 10th channel of the MCA so allowing the user to choose the window for his data collection. Section 2.2 for further details. The programme then returns return to the start

'DO YOU WISH TO RETURN TO MAIN PROGRAMME'

Exit by typing 'Y'.

#### 1.5.2) R - READ ION CHAMBERS

This command is very similar to that in the transmission menu except for the addition of two more options. When 'R' is typed at the console the VDU displays

TIME 1000

?

This allows the user to count ion chamber and fluorescence detectors for 1 second (1000 milliseconds) or enter an alteration time if desired.

The programme continues by asking the user

'OVER WHICH CHANNELS DO YOU WISH TO INTEGRATE'

These channels correspond to the multichannel analyser described above (see also 2.2).

'START=xxx'

?

The programme prints the starting channel already in memory which can be altered by inputting a new value and then proceeds to print

END=xxx

and this value may also be changed by typing in a new number. The programme then counts for the time specified at the beginning and displays the readings from both ion chambers and also the total number of counts from the multichannel analyser over the window specified.

2.1 SETTING UP AND OPERATION OF ORDER SORTING MONOCHROMATOR

Harmonic rejection is achieved with the two crystal monochromator by off-setting the first crystal out of parallel with respect to the second by a desired amount and employing a simple servo system to hold it there as the scan proceeds. Two mechanisms are provided to do this - a coarse adjustment and a fine adjustment. A schematic describing the mechanism of the order sorting monochromator is shown in fig.6. For details of the monochromator and the principles of x-ray harmonic rejection see 'An Order Sorting Monochromator for Synchrotron Radiation' G.N. Greaves, G.P. Dlakun, P.D. Quinn, M. Hart and D.P. Siddons, Nucl. Instrum. Meth. 208, (1983) 335-339 (DL/SCI/P340E).

## 2.1.1) COARSE ADJUSTMENT

The first crystal is mounted on a spring loaded lever which is moved up and down with a small stepper motor. The motor is driven by computer software by typing 'O' from the list of commands in the setting up programme and then entering the desired number of steps you wish to move the motor. The coarse adjustment is used to set the monochromator near the rocking curve peak or the required offset position from the peak. One step is approximately 0.9 arc seconds movement. Allow for two or three steps backlash when changing the motor direction. A negative number of steps decrease  $\theta$  for the first crystal.

## 2.1.2) FINE ADJUSTMENT

A small ferrite attached to the bottom of the first crystal is attracted or repelled by a solenoid fixed to the monochromator frame. Very fine adjustment of the crystal angle can therefore be achieved simply by varying the solenoid current. A 20 mA change in solenoid current is approximately equivalent to 1 step of the coarse adjust motor. The fine adjustment is controlled from the solenoid drive module located in the NIM BIN. The layout of the NIM module is shown in fig.7.

## 2.1.3) SOLENOID DRIVER

A double width NIM module is provided for driving the solenoid. With manual control selected a front panel potentiometer allows the solenoid

current to be set at any value within the range plus or minus 500 mA. The current is displayed on the left-hand DVM of the solenoid control panel. The output in millivolts equals the solenoid current in milliamps (see fig.7)

With the servo control selected, the module compares the reading from the I(O) ion chamber amplifier with a reference displayed on the right-hand DVM of the solenoid control panel (see fig.7). It adjusts the solenoid current automatically to keep the two equal. The servo reference can be set internally by a ten turn potentiometer (switch to 'INT' on the solenoid driver module) or externally by a DAC from the setting-up and scanning programmes (switch to 'EXT' on the solenoid driver module). The monochromator must first be preset sufficiently near the required set point with the coarse motor because the range of the solenoid is limited, particularly at long wavelengths. If for any reason the I(O) ion chamber reading moves out of range of the solenoid, the current will reach the 500 mA limit and the module will break the circuit. The reset button must be operated before the solenoid is used again. This is located at the bottom of the driver module. The servo mechanism operates on the 'negative side' of the rocking curve. This is defined as the side where a positive signal from the driver module increases the amount of light leaving the monochromator and hence the output of the I(O) ion chamber. The driver module will not control on the 'positive side' of the rocking curve. To locate the 'negative side' drive through the rocking curve using the coarse adjustment in a succession of negative steps. If the servo reference is set to match the I(O) output the servo can be switched on and the mechanism will control with only a few mA of current flowing through the solenoid. If this is not the case the solenoid current can be reduced by operating the coarse adjustment with the servo still switched on. A negative movement with command 'O' will drive the solenoid current less negative and vice versa.

## 2.1.4) SETTING THE MONOCHROMATOR FOR A SCAN AT 'CONSTANT I(O)'

Use the following procedures:

- 1) Move the monochromator to the start of the scan (command 'M') and peak the I(O) chamber reading as far as possible with the coarse adjust motor (command 'O').

- 2) Peak the I(O) chamber reading precisely with the solenoid driver on manual control of the driver module - NOTE THIS READING.
- 3) Calculate the I(O) chamber reading for the required offset, and set the servo reference to this value using the ten turn potentiometer with the reference switch set to 'INT'.
- 4) Zero the solenoid current, but leave the manual control selected.
- 5) Use the coarse adjust motor to bring the I(O) chamber reading as near as possible to the required set point. If the monochromator has only been exposed to beam for a short period leave for two or three minutes at this point and then re-adjust with the coarse motor. Use negative steps to come down the rocking curve to the set point (see section 2.1.3 above).
- 6) Switch the control module from manual to servo. The I(O) chamber reading should now move to within 1% of the reference value and the solenoid current should be close to zero. Small changes may be made to the reference value with the servo on. The solenoid current will alter to accommodate the change. If the solenoid current is too large (positive or negative), reduce this by operating the coarse adjustment, the servo remaining on (see section 2.1.3 above).
- 7) On some long EXAFS scans the drift in solenoid current from the start to finish may be large. The procedure in earlier paragraphs is designed to start the scan with the solenoid current near zero. If it is known that a large positive drift in I(O) will occur over the scan, then the coarse motor should be adjusted to produce a negative solenoid current (say, -100 to -150 mA), at the start of the scan. This increases the effective range of the solenoid over the scan.

#### 2.1.5) SETTING THE MONOCHROMATOR FOR A SCAN WITH 'CONSTANT' HARMONIC REJECTION

- 1) Move monochromator to the start of the scan (command 'M') and peak the I(O) chamber as far as possible with the coarse adjust motor (command 'O').

- 2) Peak the I(O) chamber reading precisely with the solenoid driver on manual control in the NIM BIN - NOTE THIS READING.
- 3) Move the monochromator to the end of the scan and peak the I(O) chamber as far as possible with the coarse adjust motor.
- 4) Peak the I(O) chamber reading precisely with the solenoid on manual control - NOTE THIS READING.
- 5) Move the monochromator back to the start of the scan. Remember to take the backlash out of the monochromator motor. This is achieved by adding 50 more than is required to reach the start position, and then move down 50 steps to the beginning of the scan.
- 6) Exit from moving the monochromator, this will return you to the menu of main commands. Type 'D' to set the percentage harmonic rejection for double crystal monochromator and reply to the inquiries (explanation of this command is dealt with in section 1.4.11 above). When all the prompts have been answered the programme returns to the main list of commands - then you should switch to the external reference.
- 7) Use the coarse adjust motor to bring the I(O) chamber as near as possible to the required set point which is now displayed on the servo reference. If the monochromator has only just been exposed to white radiation do not carry out any movement of the first crystal until the system has equilibrated. This will take approximately three minutes. Make sure when you sit at the set point on the rocking curve it is the negative side (see section 2.1.3 above).
- 8) Switch the control module from manual to servo. The I(O) chamber reading should now move to within 1% of the reference value. At each new monochromator position the reference value will be recalculated. This is obtained by linear interpolation between the start and end values of I(O) inputted in command 'D'. The method allows for a more constant harmonic rejection to be achieved. This is particularly important at short or long wavelengths where the intensity of the SRS changes rapidly with wavelength.

## 2.2 SETTING UP AND OPERATING OF FLUORESCENCE DETECTORS

At present five scintillation detectors are used, four look at the front face of the sample with the remaining looking at the back face. More detectors can be accommodated if the station is used with the mirror, as the reflected beam facilitates more space around the sample. Detectors are arranged on the 18 cm sphere and an optimum geometrical arrangement is described in DL/SCI/P380E.

The fluorescence detectors are connected to a summing amplifier which is connected to a fast multiple channel analyser (MCA). The MCA spectrum is integrated between start channel and end channel, for a particular excitation energy (Bragg angle or monochromator position). The start and end channels are defined using command 'R' (1.5.2). The integrated counts are printed in the sixth column on the console. The excitation spectrum is plotted [as MCA counts/reference counts] on the session, simultaneously with the absorption spectrum.

### 2.2.1) SETTING UP PROCEDURE

- 1) Check the individual gains of each photomultiplier tube (PMT) one tube at a time by disconnecting others from the summing amplifier.
- 2) Display the MCA spectrum on the Hytec VDU by using command 'A' (1.5.1).
- 3) If the gain of a scintillator + PMT needs adjustment, this can be achieved by trimming the HT using the potentiometer provided for the individual detector on the common power supply.
- 4) Check the harmonic content by observing that the MCA spectrum is symmetric. If a broad shoulder exists on high energy side (higher channel number), harmonic rejection is not sufficient. This can be increased using command 'D' (1.4.11).
- 5) The summed signal from all the detectors should be less than 500,000 counts/sec. This limit is defined by the detectors and amplifying electronics. Ensure that the MCA does not overflow for the counting

time as defined in the scan parameters (see section 1.5). This can be rectified by changing the lower and upper level discriminator settings. The changes required are extremely small, so care should be taken [if in doubt contact S.S. Hasnain or P.D. Quinn].

- 6) Define the common limits for integration using command 'R' (1.5.2). These should encompass the bulk of the MCA spectrum, removing any low energy pile up and any high counts due to scatter or incomplete harmonic rejection.
- 7) If detectors can not be matched in this way contact S.S. Hasnain or P.D. Quinn.

### 2.2.2) CHOICE OF DETECTOR FILTERS

Fluorescence radiation is accompanied by scattered radiation (excitation wavelength) for all the detectors. The scattered contribution is small (and not worth worrying about) for the in-plane 90 degree detector [see DL/SCI/P380E], for all other detectors it is preferable that a filter is used with the collimator which is provided. Filters are essential if the metal concentration is less than 20 milli Molar.

Material for a filter is chosen such that scattered radiation is preferentially absorbed compared to the fluorescence signal. For all 3d metals for instance, filter made from the (Z-1) element is optimum, while for higher Z atoms filters made from (Z-1) or (Z-2) are most appropriate. Table 1 illustrates the principle clearly, for five examples, with  $Z=23$  to 48, the filter material is chosen such that its K absorption edge is at a higher energy than the energy of sample fluorescence. Thus for  $23 < Z < 36$  filter material can be made from (Z-1) atom, for  $37 < Z < 47$  material can be from either (Z-1) or (Z-2) while for higher  $Z > 48$ , a choice can be made from (Z-1), (Z-2) or (Z-3). The optimum thickness of the filter depends on a number of factors including the concentration of the sample, matrix (the bulk of the sample), the size of the beam, the energy of the x-ray beam etc. As a guideline, the thickness should be chosen such that the scattered intensity (as measured before the edge) is reduced by a factor of ten but the fluorescence signal (i.e. the difference of above edge  $Sc+Sf$  and pre-edge  $Sc$ , as defined in DL/SCI/P380E) is not reduced to less than 50%. A typical thickness range might be 5-15  $\mu\text{m}$ .



		Table 1		K-edge	
Z	Element	E scatter (eV) or EXAFS range	E fluo (eV)	Possible filter	Absorption edge for filter (eV)
48	Cd	26720 - 27300	22300	Ag Pd Rh	25560 24357 23260
42	Mo	20000 - 20800	17400	Nb Zr	18990 18000
30	Zn	9600 - 10200	8600	Cu	8980
25	Mn	5640 - 7200	5900	Cr	5990
23	V	5450 - 6000	4940	Tl	4965

### 2.3 OPERATION WITH TOROIDAL MIRROR

#### 2.3.1) DESCRIPTION

A cylindrical mirror is available to enhance the intensity and brightness of the monochromatic beam at the sample. The wavelength range is limited in this mode to wavelengths  $\geq 1.1 \text{ \AA}$  i.e. Zn is the heaviest element for which this mode is useful in K-edge measurements while Ta is the heaviest element for L-edge measurements.

The mirror is Pt-coated quartz, 58 cm long, 5 cm wide designed for a glancing angle of 7 mrad with  $F1=11.5 \text{ M}$  and  $F2=5.8 \text{ M}$ . Ray tracing programme and experimental evaluation show that when the mirror is bent in the vertical plane it is able to intercept 3 mrad of horizontal aperture and all of the vertical aperture (0.2 mrad). Note that without the mirror 0.7 mrad of horizontal radiation (12 mm) is used, a limit imposed by the present monochromator crystal. If full vertical aperture is used the energy resolution is degraded simply as if the slit at the monochromator had been increased to 3 mm. Slits are provided in front of the mirror to define the vertical and horizontal acceptance of the mirror, resolution compatible with EXAFS and XANES can be easily be achieved. It is recommended that XANES and EXAFS are recorded separately with different mirror acceptance.

The gain in intensity at higher energies is small particularly for a resolution compatible with XANES but this improves as the energy is lowered, e.g. for Cu XANES where 2 eV resolution is required, intensity gain is approximately 2, while for Fe the gain is  $>4$  and for Ca the gain is  $>6$  - all compared to an unfocused 12 mm beam. For EXAFS applications, these gains in intensity can further be increased by roughly a factor of two.

The image produced by the mirror is close to 2:1 demagnification of the source at the sample position i.e. 3 mrad horizontal and full vertical beam is condensed into a spot at 17.3 M from the source of less than  $6 \times 0.5 \text{ mm}$  resulting in a brightness increase of approximately 60 over the unfocused beam. This is of course an important improvement particularly for dilute and/or small specimens if intensity can be traded for energy resolution.

#### 2.3.2) ADJUSTMENT

Installation of the mirror in the beam is fairly involved - particularly its optimisation for energy resolution and intensity gain. Adjustments to the mirror are therefore only to be done by the in-house group. Users requiring the mirror should let the laboratory know well in advance so that mirror installation can be properly scheduled. Horizontal and vertical slits at the mirror are provided with manual adjustments which are available to the user.

Finally the mirror is separated from the monochromator by a beryllium window in order that it operates under high vacuum conditions. The pressure in the mirror vessel is metred at the station and should normally be less than  $10^{-7}$  torr. If this is not the case the port 7 will normally close and cannot be re-opened until the high vacuum recovers. This may occur when the mirror is first placed in the beam. If this happens the station master and the beam port coordinator should be called.

### 2.4 GUIDELINES FOR EXAFS MEASUREMENTS

Facilities have been described for measuring x-ray absorption in transmission and fluorescence geometry. The following information is designed to help newcomers make the best use of their beam time.

#### 2.4.1) MODEL COMPOUNDS

The comparative nature of the EXAFS technique necessitates the use of model compounds. These should be chosen, if possible to be similar to the local chemistry anticipated for the unknown system. They can generally be chosen from the three main groupings of metallic, covalent and ionic lattices. Model compounds enable the fitting parameters in EXAFS analysis, notably the phaseshifts, to be reliably chosen (see section 3.3). Make sure the model compound structures are known with precision. This is not always true of crystalline systems. If a change in coordination number is expected for instance, choose model compounds that exhibit the two extreme local structures. The valence state of a transition metal can often be deduced from the location of a transition metal can often be deduced from the location of the absorption edge, in which case the chemical shifts from a range of model compounds of known valence state should be measured. In the case of dilute systems model compounds of similar dilution should be used wherever possible.

#### 2.4.2) SAMPLE THICKNESS

When carrying out transmission measurements the sample should have an absorbance  $u \cdot t$  of about 1-2 at just above the absorption edge, where  $u$  is the absorption coefficient of the material and  $t$  is the thickness of the sample (typical sample thickness for transmission EXAFS might be 5-10  $\mu\text{m}$ ). When preparing the sample, make sure it is of uniform thickness. If the sample has too high an absorption make it thinner or dilute it with boron nitride or solvent. If the sample is in the form of a powder grind the two materials - specimen and dilutant - together to a fine homogeneous powder. Liquid samples can be contained in cells but ensure the region intercepted by the beam is free of bubbles. For fluorescence measurements the optimum sample thickness depends on the atomic weight of the host material. For aqueous solutions a  $u \cdot t$  of around 0.5 is useful for concentrations  $>2 \text{ m Molar}$  and  $Z < 30$ . For concentrations of  $0.5 \text{ m Molar}$  and less a thickness of 1.5-2  $u \cdot t$  is required. Typically for FLEXAFS samples might be several millimetres of thickness.

#### 2.4.3) SAMPLE AREA

The unfocused monochromatic beam dimensions are approximately  $1 \text{ mm} \times 12 \text{ mm}$ . Specimens measuring  $5 \text{ mm} \times 20 \text{ mm}$  will require only a little alignment. For slightly smaller samples the beam dimensions can be reduced

but at the expense of x-ray intensity. For extremely small specimens a focused beam should be used (see section 2.3). This measures approximately  $0.5 \text{ mm} \times 6 \text{ mm}$  and offers a gain in x-ray intensity at longer wavelengths that can be traded for a smaller beam size still.

From the sample thickness and area the quality of material required can be judged. This is always less than 1 gm and may be as little as a few mgm.

#### 2.4.4) SAMPLE ALIGNMENT

A badly aligned sample can result in:

- 1) reduced signal in the final ion chamber because of the beam catching the sample jig;
- 2) saw-tooth spectrum where the exit beam moves periodically on and off the sample as the table height is altered;
- 3) increased signal in the final ion chamber because part of the beam is missing the sample altogether - this will result in a diminished absorption edge.

Aligning the sample can be carried out by placing a polaroid at the sample stage and marking its position in the stand. Open the shutter for approximately two seconds. Remove the polaroid and develop it. Replace the developed photograph in its original position at the sample stage. Remove the signal ion chamber. Switch on the laser and align the spot onto the centre of the beam on the photograph. Once this has been accomplished remove the polaroid and check that both ion chambers are at the correct height. Finally, align the sample with the laser.

#### 2.4.5) MOVING TO A NEW EDGE

When moving the monochromator through 10 degrees or more (command 'M'), the monochromator should be moved in stages of 5000 steps. If the table is synchronised to the movement of the beam it will move in sympathy. However, this movement is not precise over the whole monochromator range of 0 to 90 degrees. It should therefore be checked periodically by moving the table up and down (sub commands 'P' and 'G' of main command 'E') and noting if there is any change in the values of the ion chambers - a 200 step movement should be sufficient. If there is, move the table in the direction where  $I(0)$  or  $I(T)$  increases until the change disappears. Remember

that in moving to longer wavelengths the sample will become more opaque and I(T) may go to zero if the sample is left in place. Slit heights should be typically set at half a millimetre for 1 Å and two millimetres for 3 Å.

#### 2.4.6) ENTRANCE AND EXIT SLITS

The height of the entrance slits define the energy resolution of the monochromator,  $dE/E$ . This is given by

$$\frac{dE}{E} = \cot \sigma \frac{S}{P}$$

where  $\sigma$  is the Bragg angle, S is the entrance slit height and P the source to monochromator distances. For the various EXAFS stations the values of P are 7.1 (16 m), 7.4 (55 m) and 9.4 (20 m). Typically at 7.1 S should be set to around 0.5 mm for 1 Å work, 1 mm for 2 Å work and 2 mm for 3 Å work. Values of S for other stations can be found by scaling by P.

Exit slits are designed to pick up any scatter from the entrance slits that has traversed the monochromator. Except for particularly exacting experiments, exit slits should be left several millimetres wider than the entrance slits. If they are set closer than this they may catch the exit beam and some of the symptoms of poor sample alignment will be experienced.

#### 2.4.7) AMPLIFIER GAINS

Typical settings for Keithley amplifiers are in the range of  $10^8$  to  $10^{10}$  for both reference and signal ion chambers. They do not need to have the same gains, but the same rise time should be used of 300 milliseconds and there should be zero suppression. Saturation on any particular range setting results in a flat output of approximately 11 volts. Should this occur the gain should be reduced or the x-ray intensity attenuated.

#### 2.4.8) ION CHAMBER SETTINGS

The Ortec power supply which provides the HT for the ion chambers should be set at 300 V. Before running a scan make sure you have readings from both ion chambers without any beam (use command 'R'). Usual dark current readings are between 500-2000 counts in one second. If no counts or too many counts are registered, use the zero adjust on the front of the Keithley which can be altered using a small screwdriver.

#### 2.4.9) CHOICE OF FLUORESCENCE FILTERS

A variety of filters are available for covering the PMTs of the fluorescence detector array. These are mainly 3d metal foils (see section 2.2.2). Other filters should be provided by the user following the rules laid out in DL/SCI/P380E.

A typical HT setting for the fluorescence PMTs is about 1.5 KV. Before beginning a set of measurements at a particular absorption edge each PMT should be checked in turn as described in 2.2.1.

#### 2.4.10) THE DON'TS

- 1) Don't run spectra with the monochromator and slit vesael's let up to air.
- 2) Don't empty the ion chambers prior to filling without first switching off the HT.
- 3) Don't switch on the PMT HT directly without first ramping from zero. (Be gentle).
- 4) Don't drive the coarse motor of the order-sorting monochromator more than 20 steps without taking a record subsequent alignment can be tedious.
- 5) Don't alter the height of the front of the experiment (commands 'E' and 'F') without altering the height of the back of the experiment (commands 'E' and 'G').

The following description provides a guide through the EXAFS data reduction and analysis programmes. A number of programmes exist in the SRS Programme Library and here we briefly describe three of them, namely:

EXCALIB:- starting from the experimental data file, this produces normalised spectra of absorption/fluorescence versus electron energy in eV and Hartrees.

EXBACK:- takes absorption from EXCALIB and produces normalised EXAFS spectra and Fourier transforms.

EXCURVE:- curved wave and Fourier transform analysis programmes.

All three programmes are interactive on TSO and the user should be in a region equal to 650 K and using a Tektronix or Tektronix compatible terminal. The description provided is only meant as an introduction to their use - full documentation exists in the SRS Programme Library. For details from TSO, type

SRS HELP

To obtain hard copy documentation for these programmes type

RECORD

followed by

SRS EXCALIB/EXBACK/EXCURVE

Wait for details to be completed on the VDU, then type

EXIT

For fuller details on EXCURVE type

EXCURVE PRINT.

### 3.1 EXCALIB

This interactive programme reads the data as produced by the data acquisition programme and reduces it to an x-ray absorption spectrum like that shown in fig.2 or the spectra shown in fig.5. The programme is executable from a registered SRS ID by typing:

EXCALIB 'carriage return'

(In the following descriptions the programme responses are preceded by an asterisk.)

* TEKTRONIX TERMINAL? Y	This programme can also be executed from a normal VDU in which case type N.
* ENTER DS NAME 'EXP1.Rxxxxx.DATA'	
* REF. OFFSET? 150	This is the dark current reading of the ion chamber amplifier.
* TRAN. OFFSET? 250	
* GAIN RATIO? 1	This is the ratio of I(t)/I(0) amplifier gains.
* PLOT SPECTRUM? N	A response of Y is also possible here.
* READ EXAFS SPECTRA.? N	If more than one scan of a sample is to be averaged type Y, in which case the DS name for the subsequent scan will be requested.
* READ INST. FUNCTION? N	If answer Y the DS of the instrument function will be asked for. This is the spectrum obtained without the sample in place.
* MONOCHROMATOR TYPE, SI(111) (0), SI(220) (1), GE(111) (2) ? 0	
>>> CORRECTED SPECTRUM	

```

* PLOT SPECTRUM ?
Y

* DEGLITCH                               Deglitching enables instrumental ar-
Y                                           tifacts to be removed and replaced
                                           by straight lines or curves.

* KEEP OUTPUT SPECTRUM ?
Y

* ENTER DATA SET NAME
Rxxxx.DATA

* RESTART ?
N

READY

```

### 3.2 EXBACK

This interactive programme is available for removing background from either absorption or fluorescence excitation spectra. An absorption spectra as recorded may look like the example of the Ge K-edge shown in fig.8.

Pre-edge background (1) is subtracted by defining a polynomial (of up to order 3) by choosing two points, N1 and N2, in the pre-edge region. A post-edge background (2) is subtracted by defining a polynomial (of up to order 4) by choosing two points, N3 and N4. The choice of N3 is quite critical and in most cases a point chosen along the rising part of the first EXAFS peak gives a satisfactory background correction. The normalised EXAFS spectrum from fig.8 is shown in fig.9.

The programme can be executed at a Tektronix (or equivalent) terminal by typing

```

EXBACK SJ R42624A 'carriage return'
(using an SRS ID of SJ and a data set R4624A.DATA purely as an example).

```

You then proceed as follows:

```

* ENTER NUMBER OF DATA RECORDS TO BE SKIPPED
2

* ENTER COLUMN NUMBER FOR ABSCISSA VALUES
1
If the absorption spectrum was produced from EXCALIB, column 1 is energy in eV.

* ENTER COLUMN NUMBER FOR ORDINATE VALUES
2

* ENTER POINT FREQUENCY
2
This is the plotting frequency; all the points are saved in the background subtracted data. If there are more than 1024 points, use a point frequency >2.

* ENTER WEIGHTING FACTOR 0, 1, 2, 3
0

* ENTER INPUT ENERGY FLAG, 1 FOR eV,
2 FOR HARTREES
at this stage the spectrum is plotted.

* SELECT EZERO
Cross-wires will now appear on the screen. Bring the cross-wires to the correct position and press any key on the console except for '0' or 'RETURN'. It is important to get ED right as this calibrates the EXAFS spectrum.

* ENTER 1 if POINT IS OK 0 TO SELECT AGAIN
1

```

\* SELECT N1 AND N2

The cross-wires will appear again. With them choose N1 and N2 respectively. The fit is then plotted.

\* ENTER 1 IF FIT IS OK 2 TO SELECT AGAIN

1

\* ENTER 0 TO REPLOT 1 TO GO ON

1

\* SELECT N3 AND N4

Again, the cross-wires appear and you choose N3 and N4. The fit is plotted and this should approximate to the atomic absorption after the edge.

\* ENTER 1 IF FIT IS OK 2 TO SELECT AGAIN

1

\* ENTER 1 TO REPLOT 0 TO GO ON

1

\* 1 TO SELECT N5, 2 N5=EZERO

2

N5 can be chosen anywhere as the calibration has already been made. N5 is the point from which the background subtracted data will be saved and plotted. The background subtracted spectrum is then plotted from N5.

\* 1 TO GO ON, 2 CHOOSE N5, 3 CHOOSE N1,  
4 CHOOSE N3 ?

If a user is not satisfied with the background subtracted he can choose to go back to pre- or post- edge subtraction or to simply re-define N5.

READY

### 3.3 EXCURVE

This is a comprehensive interactive programme starting from a normal spectrum (e.g. fig.9). It involves curve fitting to the normalised spectrum or its Fourier transform using calculated phaseshifts and the curved wave approximation. It can be executed simply by typing:

EXCURVE 'carriage return'

The following notes demonstrate only the basic options of EXCURVE. Other aspects of this comprehensive programme can be obtained from the documentation. The data sets used below serve only as examples of the required input to the programme.

\* UNITS?

1 EV+ANGSTROMS

2 BOHR RADII+HARTREES

1

\* POINT FREQUENCY ?

1

\* DSN FOR EXPERIMENT ?

A708.EXBACK

The .DATA qualifier is not required.

\* COLUMN COMBINATION ?

32

\* NO. OF HEADER LINES ?

10

\* NO. OF DATA POINTS READ = xxx

This should be less than 200 if the iteration facility will be required later on. If not, either reduce the range or increase the point frequency.

\* DSN FOR PARAMETERS ?

'EXOUTA'

A response from you of 'w' will bring the default DS into action. Again, the .DATA qualifier is not required.

\* ENTER NUMBER OF PHASESHIFT FILES

2

The phaseshift files are usually in partitioned data sets as in this example. Here we have data set PHCA06.DATA with members SC and O.

\* ENTER DSN FOR CENTRAL ATOM

PHCA06.DATA(SC)

\* ENTER DSN FOR SECOND ATOM

PHCA06.DATA(O)

\* ENTER COMMAND

L

Command L lists the parameters. Note - all parameters are either in atomic units or in Angstroms and eV, including VPI. They come in the following tabular form:

SHELL	NI = x.xxx	TI = x.xxx	RI = x.xxx	AI = x.xxx
"	"	"	"	"
"	"	"	"	"
"	"	"	"	"

\* ENTER COMMAND

P

Command P calculates and plots the EXAFS spectrum with the experiment.

Here are some more examples of EXCURVE commands:

CP Calculates and plots EXAFS and Fourier transforms with the experiment. Note - FI, the fit index, is printed on the graph.

GS

This command sets various options, e.g.:

Device options - D 1 Tektronix, D 2 Tektronix and Versatec, D 3 Versatec.

By typing D 2 spectra will be displayed on the Tektronix and stored in a file for the Versatec.

x-axis options - X 1 energy, X 2 wavenumber.

Typing X 2 will result in plots against k.

Weighting options - W 1 none, W 2 k, W 3 k<sup>2</sup>, W 4 k<sup>3</sup>.

Typing W 2 will produce a k-weighted spectra.

Other options are available and information on these is obtainable by typing

GS

You then get the response

ENTER KEYWORD + OPTION NUMBER, LIST, RESET OR EXIT

Typing 'L' will list all the GS command options.

After making your choice or choices type E to exit out of the GS command.

IT DR This command iterates all the shell radii values (Rx) although you may choose only to refine a few distances. You will be asked

HOW MANY STEPS ?

to which an adequate reply would be 20. During the iteration process the current spectrum can be plotted on the Tektronix.

Also the spectrum can be stored in the data set PLOTFL.DATA for Versatec if device 2 is in operation (see above notes on command GS).

END

On executing this command, the programme asks

GRAPHICS DATA IN PLOTFL.DATA

VERSATEC PLOTTER

DO YOU WISH TO SUBMIT PLOTS ?

If the answer is 'yes' then

JOB SJB PLOT SUBMITTED

READY

Should you need more information whilst executing the programme, type H.

APPENDIX

Details of edge positions, spectral ranges and imitation chamber gas partial pressures for K and L<sub>III</sub> edges in the wavelength range 3.07 to 0.374 Å. Monochromator positions are given in millidegrees for Si(111) and Si(220) crystals.

5	THE VALUES QUOTED ARE FOR THE K EDGE OF THESE ATOMS				SCAN RANGE (IN TERMS OF MOTOR POSITION) FOR A TYPICAL SCAN				GAS PARTIAL PRESSURES			
10	ELEMENT	EDGE POSITION	LAMBDA IS IN ANGSTROMS		SI(220) & SI(111) ARE IN MOTOR POSITION		SI(110) SI(220) SI(111)		20%		80%	
20 :			SI(220)	SI(111)	SI(220)	SI(111)	SI(110)	SI(220)	SI(111)			
30 :			SI(220)	SI(111)	SI(220)	SI(111)	SI(220)	SI(111)	SI(220)	SI(111)	SI(220)	SI(111)
40 :			SI(220)	SI(111)	SI(220)	SI(111)	SI(220)	SI(111)	SI(220)	SI(111)	SI(220)	SI(111)
50 :	LAMBDA	KEY	SI(220)	SI(110)	SI(220)	SI(111)	SI(220)	SI(111)	SI(220)	SI(111)	SI(220)	SI(111)
60 K	3.43645	3.689	63497	33235	72310-47998	35696-27454					5.4 Torr Å	39 Torr Å
70 Ca	3.07016	4.038	53084	29318	59742-44119	31938-25236					7.1 Torr Å	51.4 Torr Å
80 Sc	2.75720	4.697	45891	26088	47249-37556	26726-21920					9.1 Torr Å	66 Torr Å
90 Ti	2.49730	4.965	40567	23472	43790-35464	25080-20814					12.1 Torr Å	87.3 Torr Å
100 V	2.26902	5.464	36220	21216	38700-32171	22515-19032					15.1 Torr Å	108.9 Torr Å
110 Cr	2.07012	5.989	32622	19279	34579-29342	20340-17464					19.2 Torr Å	138.5 Torr Å
120 Mn	1.89636	6.538	29594	17605	31171-26893	18481-16083					24.1 Torr Å	174.2 Torr Å
130 Fe	1.74334	7.112	27000	16144	28293-24750	16875-14857					30.8 Torr Å	222.5 Torr Å
140 Co	1.60811	7.710	24758	14861	25832-22864	15477-13766					38.4 Torr Å	277.1 Torr Å
150 Ni	1.48802	8.332	22799	13729	23702-21191	14252-12790					47 Torr Å	338.6 Torr Å
160 Cu	1.38043	8.982	21069	12719	21832-19692	13165-11910					56.3 Torr Å	406.4 Torr Å
170 Zn	1.28330	9.661	19523	11810	20177-18340	12195-11111					69.7 Torr Å	502.9 Torr Å
180 Ga	1.19567	10.369	18142	10993	18703-17119	11326-10386					88 Torr Å	635 Torr Å
190 Ge	1.11652	11.105	16903	10258	17387-16012	10845-9726					102 Torr Å	736 Torr Å
200 As	1.04497	11.865	15791	9594	16211-15012	9845-9128					124.3 Torr Å	97.4 Torr Å
210 Se	0.97978	12.654	14783	8990	15150-14099	9211-8580					147 Torr Å	115.5 Torr Å
220 Br	0.91994	13.477	13861	8437	14184-13260	8631-8075					179 Torr Å	140 Torr Å
230 Kr	0.86546	14.326	13025	7934	13309-12493	8105-7613					211 Torr Å	161 Torr Å
240 Rb	0.81549	15.204	12261	7473	12512-11788	7624-7188					249 Torr Å	150 Torr Å
250 Sr	0.76969	16.108	11563	7051	11785-11142	7186-6797					295 Torr Å	174 Torr Å
260 Y	0.72762	17.040	10922	6664	11051-10547	6784-6436					341 Torr Å	198 Torr Å
270 Zr	0.68877	18.001	10333	6307	10510-9996	6414-6102					406 Torr Å	231 Torr Å
280 Nb	0.65291	18.989	9789	5977	9948-9487	6074-5507					463 Torr Å	273 Torr Å
290 Mo	0.61977	20.005	9288	5673	9431-9015	5759-5507					589 Torr Å	297.5 Torr Å
300 Tc	0.58888	21.054	8821	5389	8950-8575	5467-5239					621 Torr Å	341.6 Torr Å
310 Ru	0.56047	22.121	8393	5128	8509-8170	5199-4993					725 Torr Å	386 Torr Å
320 Rh	0.53378	23.228	7990	4884	8095-7788	4948-4760					60.9 Torr Å	439 Torr Å
330 Pd	0.50915	24.351	7619	4592	7715-7435	4716-4546					70.4 Torr Å	503 Torr Å
340 Ag	0.48582	25.521	7268	4444	7355-7100	4497-4342					78.5 Torr Å	564 Torr Å
350 Cd	0.46408	26.716	6941	4245	7021-6788	4293-4151					90.3 Torr Å	656 Torr Å
360 In	0.44397	27.926	6639	4060	6712-6499	4105-3975					99.6 Torr Å	718 Torr Å
370 Sn	0.42468	29.195	6350	3884	6416-6221	3924-3805					118.4 Torr Å	639 Torr Å
380 Sb	0.40663	30.491	6079	3718	6074-5961	3755-3647					186.3 Torr Å	681 Torr Å
390 Te	0.38972	31.814	5825	3564	5881-5717	3598-3498					140.3 Torr Å	757 Torr Å
400 I	0.37379	33.169	5586	3418	5637-5487	3449-3357						

410 END



## 410 THE VALUES QUOTED ARE FOR THE L(111) EDGE OF THESE ATOMS

420 ELEMENT	EDGE POSITION	SI(220)	SI(111)	SCAN RANGE (IN TERMS OF MOTOR POSITION) FOR A TYPICAL SCAN	GAS PARTIAL PRESSURES	
430 :	LAMDA IS IN ANGSTROMS	SI(220)	SI(111)	300EV BELOW AND 600EV ABOVE THE EDGE		
440 :	SI(220) & SI(111) ARE IN					
450 :	MOTOR POSITION					
460 :	LAMDA	KEY	SI(220)	SI(111)	20%	80%
470 Mo	4.90930	2.525	51534	62714-39255		
480 Ic	4.62537	2.581	47536	56150-37063		
490 Ru	4.36632	2.840	44138	51124-35088		
500 Rh	4.13889	2.996	41308	47178-33360		
510 Pd	3.89688	3.182	38427	43324-31524		
520 Ag	3.67283	3.376	35858	40005-29824		
530 Cd	3.50070	3.542	33940	84818-51216		
540 In	3.32154	3.733	31989	70136-48172	4.7 Torr A	33.7 Torr A
550 Sn	3.16952	3.912	30365	63367-45692	5.2 Torr A	37.7 Torr A
560 Sb	2.99637	4.138	28548	57273-42958	5.8 Torr A	41.9 Torr A
570 Te	2.85162	4.384	27052	52903-40733	6.3 Torr A	45.2 Torr A
580 I	2.71901	4.560	25700	49281-38736	7.3 Torr A	53 Torr A
590 Xe	2.59330	4.781	24431	46099-36872	8.3 Torr A	60 Torr A
600 Cs	2.47227	5.015	23222	43219-35101	9.4 Torr A	68 Torr A
610 Ba	2.36114	5.251	22122	40703-33493	10.6 Torr A	76.8 Torr A
620 La	2.25792	5.491	21107	38462-32011	12 Torr A	86.6 Torr A
630 Ce	2.16410	5.729	20191	36493-30620	13.5 Torr A	97 Torr A
640 Pr	2.07707	5.969	19346	34718-29440	15.1 Torr A	108.9 Torr A
650 Nd	1.99452	6.216	18548	33077-28275	16.7 Torr A	120.4 Torr A
660 Pm	1.91888	6.461	17821	31605-27211	19.2 Torr A	138.5 Torr A
670 Sm	1.84464	6.721	17110	30188-26169	20.6 Torr A	149 Torr A
680 Eu	1.77491	6.985	16444	28880-25193	23.2 Torr A	167 Torr A
690 Gd	1.70955	7.252	15822	27674-24280	25.8 Torr A	186 Torr A
700 Tb	1.64840	7.521	15242	26560-23427	28.4 Torr A	205 Torr A
710 Dy	1.59025	7.797	14692	25510-22613	31.3 Torr A	226 Torr A
720 Hd	1.53529	8.076	14174	24533-21848	34.9 Torr A	251 Torr A
730 Er	1.48242	8.364	13676	23603-21112	38.4 Torr A	277 Torr A
740 Tm	1.43140	8.662	13197	22713-20402	42.3 Torr A	305 Torr A
750 Yb	1.38518	8.951	12763	21915-19758	47 Torr A	338.6 Torr A
760 Lu	1.34039	9.250	12344	21147-19135	51 Torr A	367 Torr A
770 Hf	1.29570	9.569	11926	20386-18512	55.5 Torr A	400 Torr A
780 Ta	1.25427	9.885	11539	19685-17935	62 Torr A	448 Torr A
790 W	1.21529	10.202	11176	19030-17392	67.2 Torr A	485 Torr A
800 Re	1.17720	10.532	10822	18394-16860	74 Torr A	532 Torr A
810 Os	1.14143	10.862	10489	17800-16361	80.6 Torr A	581 Torr A
820 Ir	1.10599	11.210	10160	17214-15866	89.3 Torr A	644 Torr A
830 Pt	1.07306	11.554	9854	16672-15406	95.2 Torr A	686 Torr A
840 Au	1.04028	11.918	9550	16135-14947	102 Torr A	736 Torr A
850 Hg	1.00944	12.282	9265	15633-14515	117 Torr A	823 Torr A
860 Tl	0.97968	12.656	8989	15148-14097	124.3 Torr A	97.4 Torr A
870 Pb	0.95112	13.036	8725	14686-13697	136 Torr A	106 Torr A
880 Bi	0.92459	13.410	8480	14258-13324	147 Torr A	115.5 Torr A
890 Po	0.89761	13.813	8231	13824-12945	158 Torr A	123 Torr A
900 At	0.87234	14.213	7997	13419-13590	172 Torr A	133 Torr A
910 Rn	0.84845	14.613	7777	13037-12253	187 Torr A	144 Torr A
920 END						

## FIGURE CAPTIONS

- Fig.1 Schematic of the EXAFS equipment
- Fig.2 Transmission spectrum of Ni foil.
- Fig.3 XANES of V metal foil with different amounts of harmonic contamination with monochromatized beam. The percentages on the right refer to  $I/I_{\max}$ , the intensity used,  $I$ , compared to the maximum available when the two crystals are aligned,  $I_{\max}$ .
- Fig.4 Arrangement of scintillation detectors for fluorescence measurements.
- Fig.5 Fluorescence excitation spectrum of  $\text{CuSO}_4$  solution obtained in a single scan of 6 sec per energy point for a 5 m molar, 1 m molar and 250  $\mu$  molar solution.
- Fig.6 Order sorting Si 220 monochromator. M is a ferrite magnet influenced by solenoid S. C is a stepper motor driven screw and provides coarse adjustment about pivot P.
- Fig.7 Layout of control rack panel for order sorting monochromator and ion chamber currents.
- Fig.8 Absorption at the Ge K-edge of liquid  $\text{GeCl}_4$  showing the pre(1) and post (2) edge polynomials used for the background subtracting programme EXBACK.
- Fig.9 Normalized EXAFS for the  $\text{GeCl}_4$  absorption spectrum shown in fig.8.




Fig.1




Fig.2




Fig.3




Fig.4




Fig.5






Fig.6

### Servo driver



### EXAFS Servo panel



### EXAFS Current monitor




Fig.7






Fig.8




Fig9



